

NEWBORN HEARING SCREENING TRAINING CURRICULUM

COMMUNICATING WITH PARENTS/CARETAKERS

NHSTC A-ABR Troubleshooting Checklist

Issue:

- Test is taking too long to complete (> 12 to 15 minutes) There is too much noise, artifact/EEG levels are too high

Possible causes:

Myogenic and/or electrical interference

- a. Myogenic noise can occur from baby moving, e.g. sucking, crying, muscle movement, neck strain
 - ♦ Ensure baby is fed, dry and relaxed and swaddled prior to screening
 - ♦ If sucking on pacifier, remove it when testing is initiated
- b. Interference from electrical noise
 - ♦ If equipment runs on electrical power change outlet and try and use a protected outlet
 - ♦ Change room location
 - ♦ If baby is hooked up to any monitors, ask nurse if they can temporarily be
 - ♦ urned off while screening is conducted
 - ♦ If equipment uses power but can also run on battery power, unplug
 - ♦ Turn off cellphones
 - ♦ Turn off any florescent lights
- c. Other things to check
 - ♦ Is sound coming out of the transducer/s (probe, couplers, etc.)?
 - ♦ Is the probe fit secure (hand's-free fit)?
 - ♦ Are the couplers fitting secure over the baby's ears?
 - If they aren't sticking or coming off try placing the knit cap over the top of the couplers to hold them into place
 - ♦ Have any of the electrodes loosened or fallen off?
 - ♦ Are all of the cables properly connected?
 - ♦ Is there too much ambient noise in the room?
 - ♦ High noise levels near screening (talking, fans, vents, monitors)?

