

GSI Corti/Welch Allyn/New EroScan OAE Screener Quick Use Guide



To Screen - 2 Button pushes!

1. Press **DOWN ARROW** to turn unit on.
2. Press the **RIGHT** or the **LEFT ARROW** to select test ear and begin the screening test

Other Quick functions:

- OFF: The unit automatically turns itself off after 1 minute of inactivity.
- SAVE a result: The unit automatically saves the results of completed tests. Will save up to 250 results

Power:

- 3 to 4 hours to fully charge the re-chargeable battery
- 1000 tests/15 hours of testing – the screening unit will screen 1000 test or for approximately 15 hours of continuous testing before it needs re-charged
- **DO NOT** fully deplete battery for maximum battery life. Charge when 5-10% of battery life remains. Battery charge is shown by a battery symbol on the display
- You **cannot** test with the unit plugged into a power outlet while it re-charges.

PE Tubes:

- To test children with PE tubes, the AutoStart needs to be disabled. This is accomplished by first inserting the probe with appropriate ear tip into the ear canal and obtaining a proper seal. To disable AutoStart at the main menu select the ear to be tested by holding down the **<LEFT or RIGHT>** ⇒arrow keys for 3 seconds until the green 'TEST' light turns off. Once the key is released, the ERO•SCAN® will calibrate and test as before

Maintenance

- Annual calibration - Complete annually to ensure accuracy and to obtain software updates etc.
- Probe tips should be replaced when they become clogged Do not attempt to clean them - they are disposable and must be replaced when they become clogged.
- Cleaning:
 1. May be wiped clean with a damp cloth using a mild antiseptic solution (e.g. cetyl chloride).
 2. Do not put excessive pressure on the clear display window or allow any utensil to puncture the display window or keypad.
 3. Do not allow any fluid to enter the device or immerse the instrument in fluids or attempt to sterilize the instrument or any of its accessories.

